

# Vital-Cath™ 72

## Closed Suction Systems

It features isolation valve which provides complete seal to the irrigation chamber by a simple switch action.



- Isolation valve by a simple switch operation
- Ergonomic suction control switch with security lock and an easy single hand operation design
- Catheter with smooth edge of tip and side holes.



Always make sure to withdraw the suction tube completely until the black ring indicator is visible before turning the valve to “OFF” position.



Suction with simultaneous ventilation is available at **ON** position.

At **OFF** position continuous ventilation is provided for patient comfort. In the same time irrigation chamber is 100% sealed allowing better cleaning results.

## Specifications

Product Code	Catheter Length (mm)	Remark	Outer Diameter	Packaging	ISO Colour
65100A / 65100AM 67100A / 67100AM	540	Endotracheal	10FR / 3.33mm	20ea/box	●
66100A / 66100AM 68100A / 68100AM	305	Tracheostomy			
65120 / 65120M 67120 / 67120M	540	Endotracheal	12FR / 4.0mm	20ea/box	○
66120 / 66120M 68120 / 68120M	305	Tracheostomy			
65140 / 65140M 67140 / 67140M	540	Endotracheal	14FR / 4.67mm	20ea/box	●
66140 / 66140M 68140 / 68140M	305	Tracheostomy			
65160 / 65160M 67160 / 67160M	540	Endotracheal	16FR / 5.3mm	20ea/box	●
66160 / 66160M 68160 / 68160M	305	Tracheostomy			



- M = with MDI port(OPTIONAL)
- All adult size Vital-Cath™ package includes disconnect wedge and corrugated adapter.